


ATP regeneration system:

This set of reactions regenerates ATP from ADP thus keeping the ATP concentration high enough for the ATP dependent reactions of the RecA.

Aberrations:

ATP – Adenosine Tri Phosphate
ADP – Adenosine Di Phosphate
CPK - Creatine Phospho Kinase
PC – PhosphoCreatine
C – Creatine

This is an equilibrium reaction:



The PCK enzyme has a random sequential mechanism – the order of addition of substrates and release of products is random.

1) Avram lysates' system:

Cell extract + 20x Creatine kinase + 10x energy mix + H₂O = 10ul / reaction
30°C, 0 – 3 hours

20x Creatine kinase = **1mg/ml** (stock = 5mg/ml)

10x energy mix = **75mM** Phosphocreatine
10mM ATP
10mM MgCl₂

2) For add-back use in Valary's Bomb lysis buffer:

Noc lysate + 10x Regeation mix + beads + H₂O = 50ul/ reaction (30°C, 30min)

10x Regeation mix: 1M Tris.HCl pH7.2	20ul (final conc. 0.2M)
1M MgCl ₂	20ul (final conc. 0.2M)
1M Phosphocreatine (CP)	20ul (final conc. 0.2M)
10mg/ml Creatine Kinase (CPK)	10ul (final conc. 1mg/ml)
0.1M ATP	10ul (final conc. 10mM)
0.2M βME	5ul (final conc. 10mM)
H ₂ O	=> 100ul (total volume)

NOTE: Modified Bomb Lysis Buffer - (β ME instead of DTT, omit EDTA)

For 500ml:

0.5 M HEPES pH7.4	40ml (final conc. 40mM)
14.3 M βME (original bottle)	70ul (final conc. 2mM)
β-gp (glycerolphosphate)	6g (final 60mM)

Other inhibitors (1000x) added upon use.